

**TK-VETKIT**



*Portable Small  
Animal Gastroscope*

**EC-250H**




**SEESHEEN | TK-VETKIT**



**TK VETKIT CO.,LTD.**


ADD:Plant 1, No. 199, Shui'an 1st Road, Xiangzhou District, Zhuhai  
TEL:+86-755-23761871  
EMAIL:tkvet@seesheen.com  
WEB:www.seesheen.com

C-EC(250H)-202308-V1.0



## CMOS+LED Technology

Non-fiber design makes it durable and easy for maintenance.




## Features




### Multi-function Buttons


① Zoom ② Capture/Record ③ Freeze




### Water/Air Nozzle



### 4-Way Deflection



### Big Channel



### HD 720P


## Portable



### 10inch Monitor



### 4inch Monitor





### USB Water/Air Pump

## Specification

Model	Insertion Tube	Working Channel	Working Length	Water/Air Nozzle	High Definition	Bending	Rotation	View Angle
EC-250H	5.8mm	2.8mm	1500mm	Yes	1280*720	U180° D130°	L120° R120°	120°

## Small Animal Application

EC-250H with its special water/air nozzle and big working channel, ideal for small animal gastroscopy.



## Large Animal Application

With extended 1.5 meters working length and big working channel, EC-250H is ideal for equine applications, including upper and lower airway endoscopy, sinuscopy, vaginoscopy and hysteroscopy, etc.

