




FLEXIBLE DIGITAL
CYSTOSCOPE

CF Series



With fiber-free design, all-electronic technology, excellent toughness, ultra-fine outer diameter.



Specifications	CF52
Outer diameter of insertion tube	Φ5.2mm
Inner diameter of instrument channel	Φ2.6mm
Field of view	120°
Depth of view field	3mm-200mm
Range of tip Bending	(U)210° (D)210°
Working length	400mm

Specifications	VS100
Imaging resolution	1280x800
Video output	HDMI
Data storage	USB Disk, SD card(Built-in 8GB)
Inner/ Video format	JPG/ BMP/ MP4
Power	DC3.7V Li-ion Battery; 100~240V; 50/60Hz

-  +86 400-0756-365 (Service)
-  +86 (0) 756 6910150
-  301 Building 4, No. 30 Feipeng Road, GMTCM Park,
Hengqin New Area, Zhuhai City, Guangdong Province, China
-  international@vision-medt.com



www.vision-medt.com